

I have a Coloplast Prosthetic Implant

To whom it may concern:

In case of emergency, please be aware that I have a Coloplast Penile Implant that **contains metal**. The implant is **MR conditional**.

See opposite side for details.

MR Conditional

- Static Magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Maximum whole-body averaged specification absorption rate (SAR) of 3W/kg for 15 minutes of scanning produced a temperature rise of less than 1.9°C



For more information, please contact Coloplast
Surgical Support Department 800.258.3476

The Coloplast logo is a registered trademark of Coloplast A/S.
© 2022-10. All rights reserved Coloplast A/S. PM-08544

